

## Heel and Toe Supports

### Aircast AirHeel™



The AirHeel™ is specifically designed to treat plantar fasciitis, Achilles tendonitis, and heel pain. Using two interconnected aircells located under the foot arch and in back of the heel, the AirHeel applies pulsating compression with every step to help reduce swelling and discomfort, and enhance circulation. The AirHeel is made of lightweight breathable fabric and is simple to apply - just slip it on like a sock and adjust the fit with one strap.

#### Ordering information

|               |       |        |       |
|---------------|-------|--------|-------|
| <b>Part #</b> | 09AS  | 09AM   | 09AL  |
| <b>Size</b>   | Small | Medium | Large |

|               | Small      | Medium      | Large      |
|---------------|------------|-------------|------------|
| Men           | Up to 7    | 7.5 - 11.5" | 12+"       |
| Women         | Up to 8.5  | 9 - 13"     | 13.5+"     |
| Ankle Circumf | 10 - 12.5" | 12 - 14"    | 13.5 - 16" |

### Bledsoe Bunion Boot



The **Bledsoe Bunion Boot** simplifies post-op rehab for patients. In the past, a special shoe, boot, or cast was worn to accommodate the surgery and to help provide stability to the foot during the six to-eight-week recovery period. Today, with the Bledsoe Bunion Boot, the patient can be safely and comfortably weight bearing right after surgery. Most patients are thrilled to avoid the use of crutches.

#### Part # and Sizing Chart

| Shoe Size |        | Bunion Boot |             |
|-----------|--------|-------------|-------------|
| Men       | Women  | Size        | Part Number |
| 2.5-5.5   | 3.5-6  | S           | AL758003    |
| 5.5-9     | 6-9.5  | M           | AL758005    |
| 9-12.5    | 9.5-13 | L           | AL758007    |

### Bunion Splint



Flannel lined foam with brushed nylon backing for contact closure. Malleable stay and strap may be adjusted to accommodate indication. Available in two styles, universal or sized.

| Part #   | Male        | Female     | Size |
|----------|-------------|------------|------|
| 79-96133 | Right <6    | Right <7.5 | S    |
| 79-96135 | Right 6-8.5 | Right 8-10 | M    |
| 79-96137 | Right 9+    | Right 10+  | L    |
| 79-96143 | Left <6     | Left <7.5  | S    |
| 79-96145 | Left 6-8.5  | Left 8-10  | M    |
| 79-96147 | Left 9+     | Left 10+   | L    |

### Silipos Wonder Spur Soft Silicone Heel Cups



Use in dress or athletic shoes for heel spurs, plantar fasciitis, knee and back pain and Haglund heel. Duo-density, medical-grade silicone ensures support and shock absorption. Softer silicone dots cushion sensitive areas.

| Part #    | Description            | Male  | Female |
|-----------|------------------------|-------|--------|
| SIL211SM  | Silicone Heel Cups, S  | 2-4.5 | 4-5.5  |
| SIL211MED | Silicone Heel Cups, M  | 5-7.5 | 6-9.5  |
| SIL211LG  | Silicone Heel Cups, L  | 8-10  | 10-12  |
| SIL211XL  | Silicone Heel Cups, XL | 10+   | 12+    |

### USA Hallux Valgus Night Splint



| Right Prod. # | Left Product # | Size   | Men's Shoe | Women's Shoe |
|---------------|----------------|--------|------------|--------------|
| 6035-RS       | 6035-LS        | Small  | Up to 6    | Up to 8      |
| 6035-RM       | 6035-LM        | Medium | 6 - 8      | 8 - 10       |
| 6035-RL       | 6035-LL        | Large  | 9+         | 11+          |

### USA Budin Toe Splints



Aligns crooked, overlapping or flexible hammer toes with elastic bands. Promotes healing after surgery.

| Product # | Size              |
|-----------|-------------------|
| 8155      | Single Toe Splint |
| 8157      | Double Toe Splint |
| 8158      | Triple Toe Splint |

## Heel and Toe Supports

### USA Hammer Toe Crests



Relieves hammer, claw, mallet and arthritic toes. Support for bent under toes while easing pressure on toe tips. 3 per package

#### Order Information

| Deluxe Product # | Economy Product # | Size        |
|------------------|-------------------|-------------|
| 8154A-RS         | 8154B-RS          | Right Small |
| 8154A-RM         | 8154B-RM          | Right Med   |
| 8154A-RL         | 8154B-RL          | Right Lg    |
| 8154A-LS         | 8154B-LS          | Left Small  |
| 8154A-LM         | 8154B-LM          | Left Med    |
| 8154A-LL         | 8154B-LL          | Left Large  |

### USA 3-Layer Toe Separators



Soft polyfoam absorbs pressure and friction between toes. Helps to maintain proper toe alignment after surgery. 12/ package.

| Product # | Size        |
|-----------|-------------|
| 8130-S    | Small       |
| 8130-M    | Medium      |
| 8130-L    | Large       |
| 8130-XL   | Extra Large |



### USA Visco-Gel Toe Spacers



Maintains proper toe alignment. Vitamin enriched gel moisturizes as ti reduces friction and irritation between toes.

| Product # | Size  |
|-----------|-------|
| 8130-S    | Small |

## Heel and Toe Support



### Tulis GEL Heel Cups

|                    |       |
|--------------------|-------|
| Standard Regular   | 10215 |
| Standard Large     | 10216 |
| Heavy Duty Regular | 10217 |
| Heavy Duty Large   | 10218 |
| Heavy Duty Small   | 10219 |

### Tulis Standard Heel Cups

|         |       |
|---------|-------|
| Regular | 10010 |
| Large   | 10020 |

### Tulis Pro Heel Cups

|         |       |
|---------|-------|
| Small   | 10211 |
| Regular | 10212 |
| Large   | 10214 |

### Tulis Rigid Heel Cups

|         |       |
|---------|-------|
| Regular | 10100 |
| Large   | 10110 |

### Viscolas Heel Spur Cushions

| Size | Part # | Men's Shoe | Women's Shoe |
|------|--------|------------|--------------|
| A    | 4838   | 2 - 4      | 4 - 6        |
| B    | 4839   | 5 - 7      | 7 - 9        |
| C    | 4840   | 8 - 10     | 10 - 12      |
| D    | 4841   | 11 - 13    |              |

### Heel Cups and Spurs

| Size    | Part # | Men's Shoe | Women's Shoe |
|---------|--------|------------|--------------|
| Small   | U4769  | Child/Yth  | Child/Yth    |
| Regular | U4770  | 2 - 7      | 4 - 8        |
| Large   | U4771  | 8 - 13     | 9 - 12       |

### Procare Achilles Wedge

4" neoprene padded foam wedge fits into most closed heel walkers. Each section is  $\frac{3}{4}$ " in height and may be a "peeled-away" during Achilles Tendon rehab.

|               |             |
|---------------|-------------|
| <b>Part #</b> | <b>Size</b> |
| 79-95130      | U           |