

Thumb Protection

Hely & Weber CMC Controller



Places the thumb in the optimum position by simultaneously abducting the thumb and distracting (seating) the CMC joint. Two sizes minimizes inventory (no rights or lefts).

Ordering Information

| PROD # | DESCRIPTION | CIRCUMFERENCE |
|---------|----------------|---------------|
| 3800XSS | CMC Controller | 4.6-6.5" |
| 3800ML | CMC Controller | 6.5-11" |

CMC sizing wrist circ. from 4.5 - 9.5"

Hely & Weber Neoprene Thumb Orthosis



Neoprene construction provides an excellent fit of the hand and thumb. Supports the MP and CMC joints while allowing full finger motion. An excellent resting splint for thumb tendonitis and de Quervain's Tendinitis fits all. Elasticized straps for comfortable compression. Malleable stay offers excellent customization.

Ordering Information

| PROD # | DESCRIPTION | CIRCUMFERENCE |
|--------|-------------------------|---------------|
| 38123 | Neoprene Thumb Orthosis | 4.6-6.5" |

Sizing:wrist circ. from 4.5 - 9.5"

Hely & Weber Santa Barbara Thumb Splint



Constructed of Kuhl™ perforated neoprene for cool, comfortable compression. Radial, palmar, and dorsal stay pockets contain a malleable aluminum stay and a spiral stay, which can be mixed and matched for desired support. Ideal for treating Gamekeeper's/Skier's Thumb, Basal Joint Arthritis, and Thumb Tendonitis. Wrap around design for one-handed application. Fits both the right and left hand. Circumferential strap fits a wide range of thumb sizes.

Ordering Information

| PROD # | DESCRIPTION | CIRCUMFERENCE |
|--------|----------------------------|---------------|
| 3842 | Santa Barbara Thumb Splint | 4.5-9.5" |

Sizing:wrist circ. from 4.5 - 9.5"

Hely & Weber Trimable Thumb Orthosis



1/8" Nylon covered neoprene provides heat and compression. Wrap around design with built in assist strap for ease of application. Fits both the right and left hand. Taped seam on thumb can be trimmed at the IP joint.

Ordering Information

| PROD # | DESCRIPTION | CIRCUMFERENCE |
|--------|-------------------------|---------------|
| 3810 | Trimable Thumb Orthosis | 4.5-9.5" |

Hely & Weber - The Knuckle Orthosis



Provides positioning for fractures or injuries to metacarples, phalanges, MCP and IP joints. Universal size, Right and Left. Straps a infinitely adjustable for fit and comfort (pain & swelling). Distal buddy strapping eliminates rotation of the finger. Padded stays are malleable. Reversible - one splint acts as a Boxer's Fracture splint (ulnar gutter), Radial Gutter splint, and Middle and Ring Finger splints.

Ordering Information

| PROD # | DESCRIPTION | OPTION |
|--------|----------------------|--------|
| 3848 R | The Knuckle Orthosis | Right |
| 3848 L | The Knuckle Orthosis | Left |

Corflex Target Sewn Thumb



Designed to provide moderate support, warmth, compression and protection for the thumb, also available with stays for added support and immobilization. Indicated for tendonitis, and thumb strains or sprains. Manufactured from premium grade latex-free neoprene, optional plastic/metal stays.

Ordering Information

| Size | UNIV |
|----------|---------|
| W/O Stay | 88-2040 |
| W/Stay | 88-2042 |

Thumb Protection

Hely Weber PhomFit Wrist Hand Thumb



Useful in treatment of de Quervain's Tendinitis, Gamekeeper Injuries, Scaphoid/Navicular Injuries, and is used for early cast removal. Aluminum removable palmar stay allows for custom fit and adjustment. Constructed of nylon outer shell with a tricot covered foam lining for comfort and durability. Stockinette tongue and dorsal hook & loop closures offer simple application. Radial thumb stay is malleable.

Ordering Information

| PART # | DESCRIPTION | SIZE | HAND |
|--------|--------------|---------|-------|
| 447XSR | PhoemFit RHT | X-SMALL | RIGHT |
| 447SR | PhoemFit RHT | SMALL | RIGHT |
| 447MR | PhoemFit RHT | MEDIUM | RIGHT |
| 447LGR | PhoemFit RHT | LARGE | RIGHT |
| 447XLR | PhoemFit RHT | X-LARGE | RIGHT |
| 447XSL | PhoemFit LFT | X-SMALL | LEFT |
| 447SL | PhoemFit LFT | SMALL | LEFT |
| 447ML | PhoemFit LFT | MEDIUM | LEFT |
| 447LGL | PhoemFit LFT | LARGE | LEFT |
| 447XLL | PhoemFit LFT | X-LARGE | LEFT |

Size info: XS= 4.5 - 5.5", S= 5 - 6.5", M= 6 - 7.5", LG= 7 - 8.5", XL= 8 - 9.5"

Hely & Weber Titan Thumb Orthosis

455



Constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior. The adjustable lacing pull fastens with one hand. Oblique thumb stay is malleable to allow for desired thumb position. Useful in treatment of de Quervain's Tendinitis, Gamekeeper Injuries, Scaphoid/Navicular Injuries, and is used for early cast removal. Aluminum removable palmar stay allows for custom fit and adjustment. Constructed of nylon outer shell with a tricot covered foam lining for comfort and durability. Stockinette tongue and dorsal hook & loop closures offer simple application. Radial thumb stay is malleable.

Ordering Information

| PROD # | DESCRIPTION | OPTION |
|--------|---------------------------|--------|
| 455 R | Titan Thumb Orthosis | Right |
| 455 L | Titan Thumb Orthosis | Left |
| 459 R | Long Titan Thumb Orthosis | Right |
| 459 L | Long Titan Thumb Orthosis | Left |

Breg Thumbster



Circumferential, malleable thumb stay provides a unique custom fit, while the figure-8 strap enhances thumb support.

Ordering Information

| Part # | Description |
|--------|-------------|
| 10266 | Thumbster |

Hely & Weber UNO WHT®



Useful in treatment of de Quervain's tendinitis, Gamekeeper injuries, Scaphoid/Navicular injuries, and is used for early removal. Nylon outer shell with tricot covered foam lining. Radial thumb stay and Aluminum palmar stay are malleable. Removable, malleable anatomically contoured. Palmar Stay. Movable Stay Pod adjusts for sizing and acts as an anchor for the hook & loop straps. Brace won't snag clothing because there is no exposed hook. Malleable Radial

Thumb Stay. Stay pocket cover prevents migration of the palmar stay. "D" rings are offset to avoid irritation to the wrist. 3 Dorsal Loop & Lock closures are adjustable and trimmable. Stockinette tongue covers length of brace, providing comfort and easy application.

Ordering Information

| PROD # | DESCRIPTION | OPTION |
|--------|-------------|--------|
| 429R | UNO WHT | Right |
| 429L | UNO WHT | Left |

UNO WHT sizing wrist circ. from 4.5 - 9.5"

Hely & Weber Modabber Thumb



1/8" Neoprene provides compression. Brace fits left or right hand to reduce inventory. Aluminum stay and steel spiral stay are interchangeable and removable. Ideal for treating Basal Joint Arthritis, Throttle Thumb, and thumb tendonitis. Elasticized patient assist strap makes one handed application easy. Taping above stay allows trimming for freedom of the IP joint.

Ordering Information

| PROD # | DESCRIPTION | WRIST |
|---------|-------------------------|----------|
| 5803SM | Modabber Thumb Orthosis | Small |
| 5803STD | Modabber Thumb Orthosis | Standard |
| 5803XL | Modabber Thumb Orthosis | X-Large |

Wrist Circ. : Sm= 4.5-5.5", Standard= 5 - 9.5" XL= 5 - 9.5"

Thumb Protection

AliMed Freedom Thumb Spica



FREEDOM® Thumb Spica is ideal for de Quervain's syndrome and arthritis affecting the MP and CMC joints. You can make adjustments to the shell with scissors to ensure that the distal edge of the thumb post doesn't interfere with IP joint motion. Black or Beige.

Sizing: Trim with scissors for a precise fit. Small/ Medium fits most women and smaller men. Medium/ Large fits most men and larger women. Specify Left or Right.

Size and Order Information

| Part # | Description |
|--------|---|
| 5782 | Freedom Thumb Spica, Beige, Right, Small/Medium |
| 5783 | Freedom Thumb Spica, Beige, Left, Small/Medium |
| 5784 | Freedom Thumb Spica, Beige, Medium/Large |
| 5785 | Freedom Thumb Spica, Beige, Left, Medium/Large |
| 5482 | Thumb Spica, Black, Right, Small/Medium |
| 5483 | Thumb Spica, Black, Left, Small/Medium |
| 5484 | Thumb Spica, Black, Right, Medium/Large |
| 5485 | Thumb Spica, Black, Left, Medium/Large |

AliMed Freedom ThumbKeeper



Skier's thumb is a tear in one of the ligaments connecting the bones of the thumb. It's particularly common among skiers. If you suffer from Skier's Thumb, Gamekeeper's Thumb, or collateral ligament injuries then the **FREEDOM® ThumbKeeper™** is the perfect solution for you. **Sizing:** Measure width of hand at MP joints from radial aspect of index finger to ulnar aspect of little finger. Small: 2-1/4" to 2-3/4" Medium: 2-3/4" to 3-1/4" Large: 3-1/4" to 3-3/4"

Size and Order Information

| Part # | Description |
|--------|---|
| 5921 | Freedom ThumbKeeper w/D-Ring, Left, Small |
| 5920 | Freedom ThumbKeeper w/D-Ring Right Small |
| 5923 | Freedom ThumbKeeper w/D-Ring, Left, Medium |
| 5922 | Freedom ThumbKeeper w/D-Ring, Right, Medium |
| 5925 | Freedom ThumbKeeper w/D-Ring, Left, Large |
| 5924 | Freedom ThumbKeeper w/D-Ring, Right, Large |
| 5102 | Low Profile ThumbKeeper, Small, Right |
| 5103 | Low Profile ThumbKeeper, Small, Left |
| 5104 | Low Profile ThumbKeeper, Medium, Right |
| 5105 | Low Profile ThumbKeeper, Medium, Left |
| 5106 | Low Profile ThumbKeeper, Large, Right |
| 5107 | Low Profile ThumbKeeper, Large, Left |

M-Brace Thumb Spica



The **M-Brace Thumb Spica** is a neoprene-free cotton velour, thumb splint that combines the original wrist splint with a universal thumb spica splint attachment that attaches quickly for use on either the right or left hand for enhanced versatility all in one product. This thumb splint features a removable palmar, dorsal, and ulnar (pink side) stays to create the desired level of wrist support and immobilization required for any injury. The thumb support is also malleable for easy positioning of the thumb for proper healing of the thumb. Skier's thumb is a tear in one of the ligaments connecting the bones of the thumb. It's particularly common among skiers. If you suffer from Skier's Thumb, Gamekeeper's Thumb, or collateral ligament injuries then the **FREEDOM® ThumbKeeper™** is the perfect solution for you. **Sizing:** Measure width of hand at MP joints from radial aspect of index finger to ulnar aspect of little finger. Small: 2-1/4" to 2-3/4" Medium: 2-3/4" to 3-1/4" Large: 3-1/4" to 3-3/4"

Size and Order Information

| Product # | T-36R | T-36E |
|----------------------------|------------------------|-------------------------|
| Description | Description | M-Brace Wrist Splint |
| Size | Regular | Extra |
| Wrist Circumference | 5.5 - 7.9" 14-20 cm | 7.9 - 10.3" 20-26 cm |

Thumb Protection

Bledsoe Wrist Lacer w/ Thumb Spica



The **Bledsoe Wrist Lacer with Thumb Spica** utilizes malleable palmer stays to provide customized positioning and covers the extensor surface of the thumb. Durable yet supple suede with perforations provide for cool comfort.

Indications: For support and symptomatic relief of Gamekeeper's Thumb, deQuervain's Syndrome, scaphoid injury, Carpel Tunnel, and sprains and strains.

Size and Order Information

| Size | Wrist Circumference | Left | Right |
|------|---------------------------|----------|----------|
| XS | under 5 1/2" (14cm) | SA800101 | SA800201 |
| S | 5 1/2" - 6 1/2" (14-16cm) | SA800103 | SA800203 |
| M | 6 1/2" - 7" (16-18cm) | SA800105 | SA800205 |
| L | 7" - 7 1/2" (18-20cm) | SA800107 | SA800207 |
| XL | 7 1/2" - 8 1/2" (20-22cm) | SA800109 | SA800209 |

Bledsoe Trigger and Keepers Thumb Support



The **Trigger & Keeper's Thumb Support** features a wrap-around, one-hand application; flexible support from a malleable stay; and fits right and left hands and a wide range of thumbs.

Indications: Provides immobilization, support and relief for Gamekeeper's Thumb, Basal Joint Arthritis and Thumb Tendonitis.

Size and Order Information

| Size | Circumference | Product Number |
|------|-----------------------|----------------|
| S | Under 6.5" (17cm) | SA801003 |
| M | 6.5" - 7.5" (17-19cm) | SA801005 |
| L | Over 7.5" (19+ cm) | SA801007 |

Procare ThumbGUARD™



Lightweight protection and stabilization for collateral ligament injuries to the thumb. Constructed with preformed plastic outer shell and thin foam inner-liner. May be trimmed with cast scissors for anatomically correct fit. Ideal for dislocations, collateral ligament injuries (Gamekeeper) and arthritis of the thumb.

Size and Order Information

| Part # | Side | Width | Size | Part # | Side | Width | Size |
|----------|-------|-------------|------|----------|------|-------------|------|
| 79-87262 | Right | 1.75"-2.25" | XS | 79-87272 | Left | 1.75"-2.25" | XS |
| 79-87263 | Right | 2.25"-2.75" | S | 79-87273 | Left | 2.25"-2.75" | S |
| 79-87265 | Right | 2.75"-3.25" | M | 79-87275 | Left | 2.75"-3.25" | M |
| 79-87267 | Right | 3.25"-3.75" | L | 79-87277 | Left | 3.25"-3.75" | L |
| 79-87268 | Right | 3.75"-4.25" | XL | 79-87278 | Left | 3.75"-4.25" | XL |

Procare Thumb SPICA



Lightweight, perforated outer foam shell with a cool, cotton-terry liner. Elastic, circumferential, contact closure straps for proper positioning of splint. Ideal for soft tissue injury to the thumb, arthritis, tendonitis, deQuervain's Syndrome.

Size and Order Information

| Part # | Length | Size |
|----------|--------|------|
| 79-87113 | 9" | S/M |
| 79-87114 | 9" | S/M |
| 79-87117 | 9" | L/XL |
| 79-87118 | 9" | L/XL |

Measurement Info: S/M = Fits most women and smaller men, L/XL = Fits most men and larger women

Thumb Protection

Procare IMAK Smart Glove w/ Thumb®



Bead-filled wrist pad protects and cushions the underside of the wrist providing flexible support while allowing full use of the hand and thumb. Constructed with breathable cotton Lycra® with two removable stays at the thumb extension for stability and support. Removable dorsal stay provides flexible support and helps raise wrist to ergonomically desired position. Reversible to fit left or right hand. Ideal for the prevention and relief of repetitive motion injuries; treatment of mild deQuervain's, Carpal Tunnel Syndrome, arthritis, and other wrist and thumb pain. Latex free.

Size and Order Information

| Part # | Length | Size |
|----------|---------------|------|
| 79-87552 | up to 2.75" | XS |
| 79-87553 | 2.75" - 3.25" | S |
| 79-87555 | 3.25" - 3.75" | M |
| 79-87557 | 3.75" - 4.25" | L |
| 79-87558 | 4.25" + | XL |

Measurement Information: Forearm Circumference

Procare Universal Thumb-O-Prene™



Lightweight, low-profile wraparound design with two stays, rigid and spiral, for comfort and custom fit to address pathology. Neoprene construction provides warmth and compression without bulk. Universal, fits right or left hand. Ideal for sprains, deQuervain's Syndrome, tendonitis, hyperextension, and arthritis.

Size and Order Information

| Part # | Size |
|----------|------|
| 79-82700 | U |

Measurement Information: Forearm Circumference

Procare Quick-Fit® WTO



Comfortable, lightweight nylon/foam laminate with adjustable dorsal pod to accommodate various size patients. Pre-contoured aluminum palmer and thumb stays may be adjusted for proper angulation to address wrist and thumb pathologies. Elastic straps with d-ring closure provide compressive support and maximum patient comfort. Universal size,. Ideal for addressing deQuervain's syndrome, Gamekeeper's thumb, Scaphoid injuries, sprains and strains.

Size and Order Information

| Part # | Side/Size |
|----------|-----------|
| 79-82480 | R, U |
| 79-82490 | L, U |

Measurement Information: Forearm Circumference