

Walking Boots

Aircast SP Walker™



The Aircast SP Walker, a short pneumatic walking brace, is ideal for patients whose limb length or calf circumference prohibits them from using a standard height walker or for injuries that are best treated with a short walker. Each SP Walker is engineered for patient comfort while providing superior injury protection and support. The Walker's low rocker sole promotes more natural ambulation for improved mobility while the generous foot base has ample room for dressings without sacrificing comfort. The lightweight semi-rigid shell has a pre-inflated anterior aircell and two adjustable aircells located at the malleoli that can be individually inflated using the hand bulb (included) for a "total contact" fit and secure support. Because of the SP Walker's adjustability, changes in limb dimension can be accommodated throughout the healing process while still maintaining stability and immobilization.

Part # and Sizing Chart

Part #	Description	Height of Brace	Maximum Foot
01A-S	Small	7.75" (19.7cm)	4.312" (11cm)
01A-M	Medium	8.25" (21 cm)	4.687" (12cm)
01A-L	Large	9.5" (24cm)	5.25" (13cm)
01A-XL	XLarge	9.5" (24cm)	5.63" (14cm)

Part #	Description	Men Shoe	Women Shoe Size
01A-S	Small	4-7	5-8
01A-M	Medium	7-10	8-11
01A-L	Large	10-13	11-15
01A-XL	XLarge	13+	15+

Accessories: 01159 Walking Brace Extension Straps (set of 3)

Aircast XP Diabetic Walker™ System

The XP Diabetic Walker incorporates a high rocker sole for maximum off-loading of the diabetic foot and an innovative dual density insole to help eliminate pressure points. This insole, when combined with the semi-rigid shell and full aircell coverage, also helps regulate shear stress. The Walker's shell/aircell/rocker design maximizes plantar unloading and provides protection and immobilization. A Hand Bulb with Pressure Gauge is included with each System to measure individual aircell compression. This ensures regulated aircell inflation for neuropathic patients. The XP Diabetic Walker System is an ideal alternative to total contact casting. The XP Diabetic Walker can be easily removed and reapplied time after time. The XP Diabetic Walker can be secured closed with the Aircast Brace-Lok™ to help monitor patient progress. Walking Brace Extension Straps are available to increase the strap length by 4" (10 cm).



Part #	Description	Men Shoe	Women Shoe
01PD-S	Small	4-7	5-8
01PD-M	Medium	7-10	8-11
01PD-L	Large	10-13	11-15
01PD-XL	XLarge	13+	15+

Part #	Description	Ht of Brace	Max. Foot Width (at toes)	Max. Calf Circ.
01PD-S	Small	12.5"	4.5"	14"
01PD-M	Medium	13.75"	5"	16"
01PD-L	Large	14.875"	5.5"	26"
01PD-XL	XLarge	17.0"	5.63 "	26"

Accessories: 01159 Walking Brace Extension Straps (set of 3)

Walking Boots

Aircast XP Walker™



Each XP Walker features a low rocker sole for added comfort and ease of ambulation, a wider foot base ensuring plenty of room for dressings without sacrificing comfort, and relocated valves to improve ease of use and patient compliance. The XP Walker has a lightweight, durable semi-rigid shell that helps support the limb while providing full-shell protection. Lining the shell are overlapping aircells similar to those used in the Aircast patented Duplex™ Ankle Brace aircell system.

See Chart Below for Size and Order Information

Aircast FP Walker™



The FP Walker provides full-shell protection with pneumatic support. Each FP Walker has a low rocker sole for added comfort and ease of ambulation, and a wider foot base with plenty of room for dressings without sacrificing comfort. The valve placement has been relocated to improve compliance protection, and ease of inflation. The FP Walker has a lightweight, durable, semi-rigid shell that supports the limb while providing protection.

See Chart Below for Size and Order Information

Ordering Information XP and FP Walkers

Part #	Description	Men's Shoe	Women's Shoe
01F-P / 01P-P	Pediatric	Up to 4	Up to 5
01F-S / 01P-S	Small	4-7	5-8
01F-M / 01P-M	Medium	7-10	8-11
01F-L / 01P-L	Large	10-13	11-15
01F-XL / 01P-XL	XLarge	13+	15+

XP and FP Walkers Specifications

Part #	Description	Height of	Max FT Wdth	Max Calf
01F-P, 01P-P	Pediatric	12.25"	4"	14"
01F-S, 01P-S	Small	12.5"	4.5"	14"
01F-M, 01P-M	Medium	13.75"	5"	16"
01F-P, 01P-P	Large	14.875"	5.5"	26"
01F-XL, 01P-XL	XLarge	17.0"	5.63"	26"

Bledsoe Adjusta Fit Boot



The **Bledsoe AdjustaFit LC Boot** is our most economical, lightweight walker. The AdjustaFit LC is designed to enhance mobility for sprains, strains, and some types of fractures. For convenience and even faster application, Bledsoe now offers the **AdjustaFit LC Pre-assembled Boot** - ready to be fitted right out of the box! LC2 is equipped with ankle/heel pad. Optional air ankle/heel pad allows for ability to adjust ankle pressure for increased comfort and provides pneumatic ankle support (LC2 only).

Part # and Sizing Chart

Shoe Size		Leg	Boot	Part Number			
Men	Women	Length	Size	LC	LC2	LC Mid-Calf	LC2 Mid-Calf
1-2.5	1-3.5	9"-11" (23-28 cm)	XS	w/cuff AL053001	AL?53051	AL056001	AL?57001
				w/bootie AL053021	AL?53061		
2.5-5.5	3.5-6	11"-13" (28-33 cm)	S	w/cuff AL053003	AL?53053	AL056003	AL?57003
				w/bootie AL053023	AL?53063		
5.5-9	6-9.5	13"-15" (33-38 cm)	M	w/cuff AL053005	AL?53055	AL056005	AL?57005
				w/bootie AL053025	AL?53065		
9-12.5	9.5-13	15"-17" (38-43 cm)	L	w/cuff AL053007	AL?53057	AL056007	AL?57007
				w/bootie AL053027	AL?53067		
12.5-17	13+	17"-19" (43-48 cm)	XL	w/cuff AL053009	AL?53059	AL056009	AL?57009

NOTE: When ordering replace the ? with one of the following: 0=ankle/heel pad 2=air ankle/heel pad

Bledsoe Achilles Boot



The **Bledsoe Achilles Boot** saves time by eliminating casting series required to stretch the achilles tendon during rehab. Aluminum boot shell and upright system that, unlike plastic boots, can be reformed to fit virtually any leg or ankle shape. Four cast inter-locking 10° urethane wedges provide fast change of foot angle.

Shoe Size		Leg Measurement	Achilles Boot	
Men	Women	Length	Size	Part Number
2.5-5.5	3.5-6	11"-13" (28-33 cm)	S	AL?36013
5.5-9	6-9.5	13"-15" (33-38 cm)	M	AL?36015
9-12.5	9.5-13	15"-17" (38-43 cm)	L	AL?36017
12.5-17	13+	17"-19" (43-48 cm)	XL	AL?36019

Shoe size (preferred) or leg measurement taken from knee joint line to medial malleolus.

NOTE: When ordering replace the ? with one of the following: 0=ankle/heel pad 2=air ankle/heel pad

Walking Boots

Bledsoe Conformer Boot



The **Bledsoe Conformer Diabetic Boot** features a fully enclosed thick, but breathable foam cocoon attached to a unique auto-molding innersole which embeds into a specially designed pre-molded mid-sole insert in an aluminum shelled walking boot. This unique **Automold™ Dual Density** foot-bed is designed to evenly distribute body weight across the entire bottom surface of the foot, and up to one inch high on all sides of the foot. The auto-molding innersole features a skinned, cleanable smooth surface that molds itself to the foot shape.

Part # and Sizing Chart

Shoe Size		Foot	Single Cuff		Double Cuff	
Men	Women	Max Width	Conformer I		Conformer II	
			Left	Right	Left	Right
4-5.5	5-6.5	3 9/16"	AL?41104	AL?41204	AL?40104	AL?40204
6-7.5	7-8.5	3 3/4"	AL?41106	AL?41206	AL?40106	AL?40206
8-9.5	9-10.5	4 1/16"	AL?41108	AL?41208	AL?40108	AL?40208
10-11.5	11-12.5	4 5/16"	AL?41110	AL?41210	AL?40110	AL?40210
12-13.5	13-14.5	4 5/8"	AL?41112	AL?41212	AL?40112	AL?40212
14-15.5	15-16.5	4 13/16"	AL?41114	AL?41214	AL?40114	AL?40214

NOTE: When ordering replace the ? with one of the following: 0=ankle/heel pad 2=air ankle/heel pad

Bledsoe Platform Boot



The **Bledsoe Platform Boot** is designed with a flat innersole for procedures where toe dorsiflexion is not permissible. Aluminum boot shell and upright system that, unlike plastic boots, can be reformed to fit virtually any leg or ankle shape.

Part # and Sizing Chart

Shoe Size		Platform Boot	
Men	Women	Size	Part Number
1-2.5	1-3.5	XS	AL?43001
2.5-5.5	3.5-6	S	AL?43003
5.5-9	6-9.5	M	AL?43005
9-12.5	9.5-13	L	AL?43007
12.5-17	13+	XL	AL?43009

NOTE: When ordering replace the ? with one of the following:
0=ankle/heel pad 2=air ankle/heel pad

Walking Boots

Breg Vectra™ Lite Tall



The Vectra Lite Walker, an economy version of the Vectra Premium, features an open-celled foam liner that is contoured for maximum patient comfort. The Vectra Lite has an open back, with the ability to purchase a heel clip accessory for a closed back application. Vectra's Dual-Axis Rocker bottom encourages a more natural gait, even when the walker is rotated externally during use. Available in five sizes, XS – XL

Part #	Description
2841x	Vectra Lite Tall, XS-XL

See Chart Below for Size and Order Information

Breg Vectra™ Premium Tall



The Vectra Premium Walker, features a contoured liner with a moisture wicking inner layer, along with BREG's patented Airmesh® material for maximum comfort and breathability. The Vectra Premium's unique feature, the removable heel clip, which accommodates either closed or open back heel needs/applications. Vectra's Dual-Axis Rocker bottom encourages a more natural gait, even when the walker is rotated externally during use. Available in five sizes, XS – XL

Part #	Description
2840x	Vectra Premium Tall, XS-XL

See Chart Below for Size and Order Information

Order Information for Breg Vectra™ Lite Tall and Vectra Tall Walkers

SIZE #	SIZE	MEN	WOMEN
X=1	XS	< 2.5	<3.5
X=2	S	2.5-5.5	3.5 - 6
X=3	M	5.5-9	6-9.5
X=4	L	9-12.5	9.5-13
X=5	XL	12.5-17	>13

Breg Venture Air™ Tall



The Venture Air Tall features reinforced, multi-ribbed, nylon shell for maximum rigidity and optimal support, along with adjustable air chambers to increase stability and decrease pain and edema by providing compression. The struts on the Venture Air Tall are contoured for a better anatomical fit and integrated with full circumferential.

Part #	Description
1212X	Venture Air, Tall, S-L
7450X	Venture Air, Tall, Liner Kit, S-L

See Chart Below for Size Information

Breg Vectra™ Lite Short



The Vectra Lite Walker, an economy version of the Vectra Premium, features an open-celled foam liner that is contoured for maximum patient comfort. The Vectra Lite has an open back, with the ability to purchase a heel clip accessory for a closed back application. Vectra's Dual-Axis Rocker bottom encourages a more natural gait, even when the walker is rotated externally during use. Available in five sizes, XS – XL.

Part #	Description
2843x	Vectra Lite Short

See Chart Below for Size Information

Ordering Information for Breg Venture Air Tall and Vectra Lite Short Walkers

SIZE #	SIZE	MEN	WOMEN
X=1	XS	< 2.5	<3.5
X=2	S	2.5-5.5	3.5 - 6
X=3	M	5.5-9	6-9.5
X=4	L	9-12.5	9.5-13
X=5	XL	12.5-17	>13

Walking Boots

Breg Vectra™ Premium Short



The Vectra Premium Walker, features a contoured liner with a moisture wicking inner layer, along with BREG's patented Airmesh® material. It has a removable heel clip accommodates either closed or open back heel needs/applications. Vectra's Dual-Axis Rocker bottom encourages a more natural gait, even when the walker is rotated externally during use. Available in five sizes, XS – XL.

Part #	Description
2842x	Vectra Premium Short XS-XL

See Chart Below for Size and Order Information

Breg Venture™ Air Short



The Venture Air Short features reinforced, multi-ribbed, nylon shell for maximum rigidity and optimal support. The struts are contoured for a better anatomical fit and integrated with full circumferential straps for stability and support. Adjustable air chambers increase stability and decrease edema. Ideal for stable foot fractures and other foot injuries.

Part #	Description
1213X	Venture Air, Short, S-L
7011X	Venture Air, Short, Liner Kit, S-L

See Chart Below for Size and Order Information

Ordering Information for Breg Premium Short and Venture Air Tall Walkers

SIZE #	SIZE	MEN	WOMEN
X=1	XS	< 2.5	<3.5
X=2	S	2.5-5.5	3.5 - 6
X=3	M	5.5-9	6-9.5
X=4	L	9-12.5	9.5-13
X=5	XL	12.5-17	>13

Breg Controlled Range of Motion Walker



The Controlled Range of Motion Walker (CROM) is adjustable range of motion walker that offers different varieties of dorsiflexion and plantarflexion. The rocker bottom makes for near normal gait.

Part #	Description
0930X	Controlled Range of Motion Walker (CROM), S-L
70123	CROM, Liner Kit, S/M
70124	CROM, Liner Kit, L

See Chart Below for Size and Order Information

Breg Pin Cam Walker



The Controlled Range of Motion Walker (CROM) is adjustable range of motion walker that offers different varieties of dorsiflexion and plantarflexion.

Part #	Description
0955X	Pin Cam Walker, S-L
70123	Pin Cam Walker, Liner Kit, S/M
70124	Pin Cam Walker, Liner Kit, L

See Chart Below for Size and Order Information

Order Information for Breg Controlled Range of Motion and Pin Cam Walkers

SIZE #	SIZE	MEN	WOMEN
X=2	S	<6	<7.5
X=3	M	6.5-11	8-12.5
X=4	L	11.5-15	13+

Walking Boots

Procare MaxTrax™ Walker



Designed with a wider foot bed for comfort and stability following lower leg, foot and ankle trauma and post-op procedures. A low profile, rocker bottom helps promote natural gait and reduce plantar pressures. Cushioned inner and outer sole provides shock absorption to help improve patient comfort during ambulation. Soft nylon/foam liner ensures no seams are in critical post-surgical areas. Malleable uprights breakdown on standard height walkers to accommodate shorter tibia/fibula. Ideal for metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable foot fractures.

Part #	Shoe Size Female	Shoe Size Male	Size
79-95322	≤4	N/A	PED/XS
79-95323	4.5-6	≤5	S
79-95325	6.5-11	5.5-10	M
79-95327	11.5-14.5	10.5-13.5	L
79-95328	N/A	14-17	XL
79-95322	≤5 JR. Size		Ped/XS

Accessories		
79-95182	Replacement Liner	PED/XS
79-95183	Replacement Liner	S
79-95185	Replacement Liner	M
79-95187	Replacement Liner	L
79-95188	Replacement Liner	XL

Procare MaxTrax™ Diabetic Walker



Designed with a wider foot bed for added stability and a soft nylon/foam liner with no seams in critical post-surgical areas. Tri-laminate insole is designed with Impax Grid technology to redistribute load away from ulcers with removal of Impax foam blocks. Low profile, rocker bottom Cushioned outer sole for shock absorption. Optional Toe Cover and a set of Upright Pads for greater patient comfort. Ideal for replacement of total contact casting; treatment of pre-ulcerative or ulcerative conditions of the plantar surface of the foot; charcoat foot.

Part #	Male	Female	Size
79-95453	N/A	4.5-6	S
79-95455	5.5-10	6.5-11	M
79-95457	10.5-13.5	11.5-14.5	L
79-95458	14-17	N/A	XL

Part #	Part #	Gender	Size
79-95463	Impax Foam Insole		S
79-96465	Impax Foam Insole		M
79-95467	Impax Foam Insole		L
79-95468	Impax Foam Insole		XL

Accessories

79-95450	MaxTrax Diabeti Toe Cover	Size U
79-95460	MaxTrax Upright Pads (1pair)	Size U

Procare Nextep™



Durable, lightweight, design with shorter one-piece wrap-around posterior shell. Provides immobilization for acute ankle injuries and stable fractures of the foot.

Part #	Shoe Sz Female	Shoe Sz Male	Size
79-95083	3.5-7.5	4.5-8.5	S
79-95085	7-10	8-11.5	M
79-95087	10-13	11-14	L

Walking Boots

Procare MaxTrax™ ROM and Air ROM Walker



Provides protected range-of-motion following trauma or post-operative procedure. Ideal following Achilles tendon repair for controlled rehab. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar- and dorsi-flexion. Malleable uprights breakdown on standard height to accommodate shorter tibia/fibula. Cushioned inner/outer sole. Low profile, rocker bottom.

Size and Order Information- MaxTrax ROM

Part #	Shoe Size Female	Shoe Size Male	Size	Accessories		
79-95332	≤4	N/A	PED/XS	79-95192	Replacement Liner	PED/XS
79-95333	4.5-6	≤5	S	79-95193	Replacement Liner	S
79-95335	6.5-11	5.5-10	M	79-95195	Replacement Liner	M
79-95337	11.5-14.5	10.5-13.5	L	79-95197	Replacement Liner	L
79-95338	N/A	14-17	XL	79-95198	Replacement Liner	XL
79-95332	≤5 JR. Size		Ped/XS			

Size and Order Information- MaxTrax AIR ROM

Part #	Shoe Size Female	Shoe Size Male	Size	Accessories		
79-95432	≤4	N/A	PED/XS	79-95362	Replacement Liner	PED/XS
79-95433	4.5-6	≤5	S	79-95363	Replacement Liner	S
79-95435	6.5-11	5.5-10	M	79-95365	Replacement Liner	M
79-95437	11.5-14.5	10.5-13.5	L	79-95367	Replacement Liner	L
79-95438	N/A	14-17	XL	79-95368	Replacement Liner	XL
79-95432	≤5 JR. Size		Ped/XS			

Procare MaxTrax™ Ankle Walker and Air Ankle



Designed with a wider foot bed for comfort and stability following lower leg, foot and ankle trauma and post-op procedures. A low profile, rocker bottom helps promote natural gait and reduce plantar pressures. Cushioned inner and outer sole provides shock absorption to help improve patient comfort during ambulation. Soft nylon/foam liner ensures no seams are in critical post-surgical areas. Malleable uprights breakdown on standard height walkers to accommodate shorter tibia/fibula. Ideal following metatarsal fractures. bunionectomies. acute ankle sprains. stress and stable foot fractures.

Size and Order Information- MaxTrax Ankle Walker

Part #	Shoe Size Female	Shoe Size Male	Size	Accessories		
79-95342	≤4	N/A	PED/XS	79-95182	Replacement Liner	Ped/XS
79-95343	4.5-6	≤5	S	79-95183	Replacement Liner	S
79-95345	6.5-11	5.5-10	M	79-95185	Replacement Liner	M
79-95347	11.5-14.5	10.5-13.5	L	79-95187	Replacement Liner	L
79-95348	N/A	14-17	XL	79-95188	Replacement Liner	X
79-95342	≤5 JR. Size		Ped/XS			

Size and Order Information- MaxTrax Air Ankle Walker

Part #	Shoe Size Female	Shoe Size Male	Size	Accessories		
				Part #	Description	Size
79-95422	≤4	n/a	Ped/XS	79-95372	Replacement Liner	Ped/XS
79-95423	4.5-6	≤5	S	79-95373	Replacement Liner	S
79-95425	6.5-11	5.5-10	M	79-95375	Replacement Liner	M
79-95427	11.5-14.5	10.5-13.5	L	79-95377	Replacement Liner	L
79-95428	n/a	14-17	XL	79-95378	Replacement Liner	XL
79-95422	≤5 Jr. Size		Ped/XS			

Walking Boots

Procare Nextep™ Contour and Contour Air



Provides ultimate protection and edema control following trauma or post-operative procedures. A medial/lateral dual air cell system inflates both sides of the liner equally and is easily adjusted with an integrated pump for proper compression, comfort, and custom fit. A secure, full wrap-around shell provides maximum support for ligament and tendon repair, osteotomy, pinning procedures, stress fractures or post-operative bunion surgery. Cushioned, non-slip rocker sole. Washable brush nylon foam liner with toe cover for maximum comfort.

Size and Order Information- Nextep Contour

Part #	Shoe Size Female	Shoe Size Male	Size	Accessories		
79-95063	4.5-8.5	3.5-7.5	S	79-95163	Replacement Liner	S
79-95065	8-11.5	7-10	M	79-95165	Replacement Liner	M
79-95067	6.5-11	5.5-10	L	79-95167	Replacement Liner	L
				79-95440	Toe Cover 3 Pk	U

Size and Order Information- Nextep Contour Air

Part #	Shoe Size Female	Shoe Size Male	Size	Accessories		
79-95273	4.5-8.5	3.5-7.5	S	79-95195	Replacement Liner	Ped/XS
79-95275	7.5-10	8.5-11	M	79-95440	Toe cover	3 pk
79-95277	11-14	10-13	L			

USA Walker- Standard



USA Walker- Ankle



The low profile design provides comfortable stability for acute ankle sprains, soft tissue and stress fractures to the lower leg and ankle, and stable fractures of the ankle. Other features include fixed, malleable uprights, a breathable fiber leg wrap, removable heel pads, and non-skid rocker sole.

Size	Part #	Men's Shoe	Women's Shoe
XSmall	12001	up to 4	up to 5
Small	12003	4 to 6	5 to 7
Medium	12005	6 to 10	7 to 11
Large	12007	10 to 12	11 to 13
XLarge	12009	12 plus	13 plus

Size	Part #	Men's Shoe	Women's Shoe
XSmall	12011	up to 4	up to 5
Small	12013	4 to 6	5 to 7
Medium	12015	6 to 10	7 to 11
Large	12017	10 to 12	11 to 13
XLarge	12019	12 plus	13 plus

USA Air Walker- Standard



USA Air Walker- Ankle



The low profile design provides comfortable stability for acute ankle sprains, soft tissue and stress fractures to the lower leg and ankle, and stable fractures of the ankle. Other features include fixed, malleable uprights, a breathable fiber leg wrap, removable heel pads, and non-skid rocker sole. Adjustable medial/ lateral air cell system inflates both sides equally to accommodate changes in edema during rehab and maximize patient compliance.

Size	Part #	Men's Shoe	Women's Shoe
XSmall	12101	up to 4	up to 5
Small	12103	4 to 6	5 to 7
Medium	12105	6 to 10	7 to 11
Large	12107	10 to 12	11 to 13
XLarge	12109	12 plus	13 plus

Size	Part #	Men's Shoe	Women's Shoe
XSmall	12111	up to 4	up to 5
Small	12113	4 to 6	5 to 7
Medium	12115	6 to 10	7 to 11
Large	12117	10 to 12	11 to 13
XLarge	12119	12 plus	13 plus

Walking Boots

USA ROM Walker



Provides protected range-of-motion with settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights will lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar and dorsi-flexion. Lightweight, low-profile design with nonskid rocker bottom sole. Breathable, fiber liner with removable heel pad for maximum comfort.

Size	Part #	Men's Shoe	Women's Shoe
Small	12203	4 to 6	5 to 7
Medium	12205	6 to 10	7 to 11
Large	12207	10 to 12	11 to 13

USA Diabetic Walker



Intended to replace the total contact casting treatment of pre-ulcerative or ulcerative conditions of the plantar surface of the foot. Other indications may include severe Plantar Fasciitis or other foot conditions that require and would benefit from even distribution of body weight on the plantar surface of the foot.

Size	Part #	Men's Shoe	Women's Shoe
Small	D12003	4 to 6	5 to 7
Medium	D12005	6 to 10	7 to 11
Large	D12007	10 to 12	11 to 13