

## Finger Splints

### USA Finger/Thumb Splint



Vinyl/foam construction with adjustable contact closure straps for proper fit. Malleable stay may be shaped to hold fingers or thumb in exact degree of flexion or extension.

#### Ordering Information

Size	Part #	Length x Width
Small	43003	7.5" x .375"
Medium	43005	9 "x .375"
Large	43007	10" x .375"

### USA Plastic Finger Splint



#### Ordering Information

Size	Part #	Diameter	Length
#1	621	.6"	2.08"
#2	622	.66"	2.04"
#3	623	.74"	2.30"
#4	624	.75"	2.475"
#5	625	.75"	2.31"
#5.5	625.5	.9"	2.4"
#6	626	.9"	2.53"
#7	627	1.0"	2.47"

### Hely & Weber Dynadigit w/ Modabber



Dynamic Finger Extension splint for one or more fingers. Neoprene generates dynamic action at any joint. Apply hook dorsally to restrict MCP motion. Interchangeable Stays fit in Dorsal Pocket. Malleable Stay immobilizes. Coiled Stay assists extension or resists flexion. Easily attaches to the Modabber neoprene wrist cuff for infinitely adjustable.

#### Ordering Information

PROD #	DESCRIPTION	WRIST
5821R	Modabber Wrist	Right
5821L	Modabber Wrist	Left
605A	Dynadigit	Universal

Modabber sizing Wrist circ. from 4.5 - 9.5"

### Hely & Weber Auerbach Mallet Splint




Brush nylon with hook & loop contact closure. Supports the DIP joint in a neutral, extended position to prevent mallet deformity. Malleable aluminum stay can be bent to hyperextend the DIP joint. Allows full motion of the PIP joint. Soft, scaled rubber of the PIP joint. Soft, scaled rubber underside prevents migration and slippage, even when wet. Easy to wash and re-apply. Sold in units of three.

#### Ordering Information

PROD #	DESCRIPTION	SIZE
600S	Auerbach Mallet Splint	Small
600M	Auerbach Mallet Splint	Medium
600L	Auerbach Mallet Splint	Large
600XL	Auerbach Mallet Splint	X-Large

Sm= 1.4", Med= 1.8", Lg= 2.2", XL- 2.5" Tip to 2nd joint

### Hely & Weber Frap Strap



Seams are Melco taped instead of sewn to eliminate patient chaffing. Rubberized backing holds the Frap Strap in position. Washable. One size fits all. Just trim to fit.

#### Ordering Information

PROD #	DESCRIPTION	OPTION
604A	Frap Strap	Universal

### Hely & Weber - The Knuckle Orthosis



Provides positioning for fractures or injuries to metacarples, phalanges, MCP and IP joints. Universal size, Right and Left. Straps a infinitely adjustable for fit and comfort (pain & swelling). Distal buddy strapping eliminates rotation of the finger. Padded stays are malleable. Reversible - one splint acts as a Boxer's Fracture splint (ulnar gutter), Radial Gutter splint, and Middle and Ring Finger splints.

#### Ordering Information

PROD #	DESCRIPTION	OPTION
3848 R	The Knuckle Orthosis	Right
3848 L	The Knuckle Orthosis	Left

## Finger Splints

### Corflex Boxer Splint



Anatomically designed to provide comfortable support, protection and immobilization of the 4th and 5th metacarpals. Indicated for treatment of acute or sub-acute fractures of the 4th and 5th metacarpals commonly referred to as "boxer fractures". Constructed of lightweight low-density polyethylene shell with soft foam liner and palmer pad and contact

Product Information				
	S	M	L	XL
Right	37-1201	37-1202	37-1203	37-1204
Left	37-1211	37-1212	37-1213	37-1214
Size Information				
	S	M	L	XL
Wrist	4"-5"	5"-6"	6"-7"	
4th & 5th MCP	2 3/4"-3 1/4"	3 1/4"-3 3/4"	3 3/4"-4 1/4"	4 1/4"-4 3/4"

### Frog Splints



Part #	Size..
79-71963	12 pk S
79-71965	12 pk M
79-71967	12 pk L

### Finger Protectors



Part #	Waist	Size.
79-71883	24-30"	S
79-71884	30-36"	M
79-71885	36-42"	L
79-71886	42-48"	XL

### Finger Cots



Part #	Size
79-71903	1-1/2" S
79-71905	2-1/4" M
79-71907	3" L
79-71908	3-3/4" XL