

Shoulder Abduction Slings

DonJoy UltraSling III AB



Both versions of the UltraSling II provide immobilization for rotator cuff repairs, capsular shifts, Bankhart repairs, glenohumeral dislocations/subluxation and soft tissue repairs.

Features and Benefits

- Breathable extra padded fabric for greater comfort.
- Encourages effective healing by allowing the shoulder and arm to remain in a neutral position.
- Helps prevent post-operative internal rotation contractures.
- Promotes axillary air exchange to reduce risk of secondary infections.

Order Information

Part Number: 11-0450-X-06000 11-1530	Description Fits right or left De-rotation strap
---	---

See Size Chart Below

DonJoy UltraSling ER



For use after initial anterior shoulder dislocation or post operatively when external rotation is desired. The anterior part of the capsule hangs loosely during internal rotation and the labrum is able to displace medially. In external, the subscapularis tendon tightens thus closing the anterior joint cavity.

Features and Benefits

- 15° and 30° versions to suit your protocol
- Comfortable foam pillow for patient compliance
- Derotational shoulder strap to keep product in place
- Quick release buckles make it easy to take on and off
- Popular

Order Information

Part #	Description	Part #	Description
11-1343-x-06000	15° Right	11-1345-x-06000	30° Right
11-1344-x-06000	15° Left	11-1346-x-06000	30° Left
11-9100	15 deg. wedge	11-9101	15 deg. wedge 15 pk.

See Size Chart Below

Size Chart for DonJoy UltraSling III AB and UltraSling ER

x=2=S	up to 11"
x=3=M	11" to 13"
x=4=L	13" to 15"
x=5=XL	over 15"

Measure elbow crease to base of index finger

Breg Neutral Wedge



BREG's Neutral Wedge positions the arm in a natural state of external rotation, while adding the necessary tension to the anterior capsule following anterior shoulder dislocations or Bankart lesions; thereby reducing the recurrence of future injuries. Two sizes (S/M and L/XL). Universal left or right. Easy to fit pillow, sling and straps. Holder for pain Management device or beverage.

Order Information

Part # 11917 11918 11920 11919	Description Neutral Wedge, S/M (<14.5") Neutral Wedge, L/XL (>14.5") Night Pillow, Neutral Wedge, S/M Night Pillow, Neutral Wedge, L/XL
---	--

Measurement taken from olecranon to end of fist.

Corflex Ultra Shoulder Abduction Pillow w/ Sling



Features extended pillow design with malleable stay for proper positioning and V-lock strapping system to prevent rolling and shifting. Indicated for treatment after Bankart repairs, rotator cuff repairs,

Order Information

Product #				
	S	M	L	XL
	23-1901	23-1902	23-1903	23-1904
Size Guide				
	S	M	L	XL
Hip	24"-30"	30"-36"	36"-42"	42"-48"
Forea	9"-12"	12"-15"	15"-18"	18"-2

Measurement taken from hip circumference and length of forearm from wrist to elbow

Bledsoe ARC Shoulder Brace



The **Bledsoe ARC (Abduction & Rotation Control)** provides comfortable abduction positioning of the affected arm as well as fine adjustment of rotation. It enables the extremity to be positioned in neutral rotation or varying degrees of internal and external rotation. The ARC eliminates the pillow solution, using elegant, simplistic design engineering for improved function and comfort.

Product Features and Benefits

One size fits all...cuff folds back to create long, medium or short sleeve. Easily adapts for right or left arm. Aluminum waist frame is easily formed for accurate positioning & comfort. Ranges from 20° internal rotation to 15° external rotation (up to 70° with ARC-XR). Allows both abduction and rotation. Exercise hand grip keeps arm from sliding forward and reduces post-operative swelling and pain. Ratchet release waist clip fine-tunes fit and comfort. Half the weight of a pillow solution. Easy enough to be

Part #	Cool Version	Product	Size
AE050100	AE055100	ARC (Abduction & Rotation Control)	One-size-fits-all
AE050200	AE055200	ARC-XR (External Rotation) <i>XR is indicated for chronic dislocations</i>	One-size-fits-all

Shoulder Abduction Slings

Breg Slingshot Neutral



Designed to place the humerus in a neutral position of rotation, BREG's SlingShot Neutral is a cool, breathable shoulder immobilizer that includes a neutral pillow and exercise ball. Comprised of **Airmesh®**, a moisture wicking fabric, for enhanced breathability and comfort. Quick release buckles for easy adjustment. Universal left or right. Three sizes.

Order Information

Part #	Description
018X0	SlingShot Neutral, S-L

[See Size Chart Below](#)

Breg Slingshot 2



The SlingShot 2 features a comfortable, breathable **Airmesh®** sling and a 15° abduction pillow. Quick release shoulder and waist strap buckles make this product easy to apply. Includes exercise ball to stimulate circulation. Comprised of **Airmesh®**, a moisture wicking fabric, for enhanced breathability and comfort. Quick release buckles for easy adjustment. Universal left or right. Three sizes.

Order Information

Part #	Description
0850X	SlingShot 2

[See Size Chart Below](#)

Size Chart for SlingShot Neutral and SlingShot 2

x=2=S	11.5" - 12"
x=3=M	13.5" to 14.5"
x=4=L	15" to 16"

Measurement taken from olecranon to knuckles

Breg Shoulder Abduction Pillow



BREG's Shoulder Abduction Pillow is a shoulder immobilizer designed for a high degree of abduction. Inflatable bladder for abduction of 10° - 75°. Limits posterior shift of the shoulder. Universally sized abduction pillow. Universal left or right.

Order Information

Part #	Description
01851	Shoulder Abduction Pillow (Universal)