# **Finger Splints**

### **USA Finger/Thumb Splint**

# Vinyl/foam construction with adjustable contact closure straps for proper fit. Malleable stay may be shaped to hold fingers or thumb in exact degree of flexion or extension. Ordering Information Size Part # Length x Width Small 43003 7.5" x .375" Medium 43005 9 "x .375" Large 43007 10" x .375"

# USA Plastic Finger Splint



### **Ordering Information**

Size	Part #	Diameter	Length
#1	621	.6"	2.08"
#2	622	.66"	2.04"
#3	623	.74"	2.30"
#4	624	.75"	2.475"
#5	625	.75"	2.31"
#5.5	625.5	.9"	2.4"
#6	626	.9"	2.53"
#7	627	1.0"	2.47"

# Hely & Weber Dynadigit w/ Modabber



Dynamic Finger Extension splint for one or more fingers. Neoprene generates dynamic action at any joint. Apply hook dorsally to restrict MCP motion. Interchangeable Stays fit in Dorsal Pocket. Malleable Stay immobilizes. Coiled Stay assists extension or resists flexion. Easily attaches to the Modabber neoprene wrist cuff for infinitely adjustable.

### **Ordering Information**

PROD #	DESCRIPTION	WRIST	
5821R	Modabber Wrist	Right	
5821L	Modabber Wrist	Left	
605A	Dynadigit	Universal	

Modabber sizing Wrist circ. from 4.5 - 9.5"

# Hely & Weber Frap Strap



Seams are Melco taped instead of sewn to eliminate patient chaffing. Rubberized backing holds the Frap Strap in position. Washable. One size fits all. Just trim to fit.

**Ordering Information** 

PROD #	DESCRIPTION	OPTION	
604A	Frap Strap	Universal	

# Hely & Weber Auerbach Mallet Splint



Brush nylon with hook & loop contact closure. Supports the DIP joint in a neutral, extended position to prevent mallet deformity. Malleable aluminum stay can be bent to hyperextend the DIP joint. Allows full motion of the PIP joint. Soft, scaled rubber of the PIP joint. Soft, scaled rubber underside prevents migration and slippage, even when wet. Easy to wash and reapply. Sold in units of three.

### **Ordering Information**

PROD #	DESCRIPTION	SIZE	
600S	Auerbach Mallet Splint	Small	
600M	Auerbach Mallet Splint	Medium	
600L	Auerbach Mallet Splint	Large	
600XL	Auerbach Mallet Splint	X-Large	

Sm= 1.4", Med= 1.8", Lg= 2.2", XL- 2.5" Tip to 2nd joint

## Hely & Weber - The Knuckle Orthosis



Provides positioning for fractures or injuries to metacarples, phalanges, MCP and IP joints. Universal size, Right and Left. Straps a infinitely adjustable for fit and comfort (pain & swelling). Distal buddy strapping eliminates rotation of the finger. Padded stays are malleable. Reversible one splint acts as a Boxer's Fracture splint (ulnar gutter), Radial Gutter splint, and Middle and Ring Finger splints.


### **Ordering Information**

PROD #	DESCRIPTION	OPTION
3848 R	The Knucke Orthosis	Right
3848 L	The Knucke Orthosis	Left

Customer Service: customerservice@kinetichealthservices.com

# **Finger Splints**

	Corflex	Boxer S	Splint		
Anatomically designed to provide comfortable support, protection and immobilization of the 4th and 5th metacarpals. Indicated for treatment of acute or sub- acute fractures of the 4th and 5th metacarpals commonly referred to as "boxer fractures". Constructed of lightweight low-density polyethylene shell with soft foam liner and palmer pad and contact					
Product Information					
	S	М	L	XL	
Right	37-1201	37-1202	37-1203	37-1204	
Left	37-1211	37-1212	37-1213	37-1214	
Size Information					
	S	М	L	XL	
Wrist	4"-5"	5"-6"	6"-7"		
4th & 5th MCP	2 3/4"-3 1/4"	3 1/4"-3 3/4"	3 3/4"-4 1/4"	4 1/4"-4 3/4'	



**Finger Protectors** 

